Instructions for Use (EN)

Introduction

The following is a general procedure for the use of Genea Biomedx Gems In Vitro Fertilisation (IVF) medium products. These products should be used as per their intended use only. They should not be used or substituted for other applications.

Each laboratory should establish its own procedures that are optimised for the individual medical facility, including suitable control measures if using open culture to prevent evaporation of media and resultant osmolality changes.

Specifications and Quality Assurance

Each lot of product (with the exception of VitSol 3 (DMSO) in the Vitrification Set in which the tests are not applicable) is tested for:

- Sterility (Sterile Filtered)
 - No growth
- Endotoxin by LAL (limulus amoebocyte lysate) test
 - Endotoxin level <0.4 EU/mL
- Biocompatibility by MEA (mouse embryo assay) test
 - 1-cell \geq 80% of 2PN's developed to blastocyst stage or beyond by 96 hours (with the exception of Sperm Wash Gradients in which MEA is not applicable)
- Osmolality (with the exception of Vitrification Set and Gavi Medium Cartridge in which it is not applicable; refer to individual Set general information below for additional details)
- pH (with the exception of VitSol3 of the Vitrification Set in which it is not applicable; refer to individual Set general information below for additional details)
- SPTT (Sperm Toxicity Test)
 - Sperm Survival Tested
 - Applicable only to all direct sperm related media (Sperm Buffer, Sperm Medium and Sperm Wash Gradients)

All results are provided on lot specific Certificate of Analysis, available upon request.

Storage and Stability

When stored as directed, Genea Biomedx Gems IVF medium products are stable until the expiration date shown on the product label. These products cannot be re-sterilised after opening.

Discard after use.

Do not use the product if:

- Packaging appears damaged or seal is broken.
- Solution appears turbid.
- Expiry date has been exceeded.

Precautions

Where Human Serum Albumin is a component of the medium:

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens.

Source material from which this product was derived was found to be negative when tested by the manufacturer for anti-HIV 1/2, anti-




HCV, HBsAg, syphilis, HCV RNA, HIV-1 RNA, HAV RNA, HBV DNA and Parvovirus B19 DNA as per the Certificate of Analysis.

There are no reports of virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time that a Genea Biomedx Gems IVF medium product is administered to a patient, the name and batch number of the product is recorded in order to maintain a link between the patient and the batch of the product.

Contraindications

Where gentamicin is a component of the medium do not use on a patient with a known allergy to gentamicin or similar products.

Definition of Symbols

	Manufacturer
	Batch Code
	Reference
	Use by
	Sterile using aseptic processing technique
	Do not re-sterilise
	Do not re-use
	Caution. Consult instructions for use
	Temperature limitation
	Keep dry
	Keep out of sunlight / Keep away from heat
	This way up
	Fragile, handle with care
	Do not use if package is damaged
	Federal Law (USA) restricts this device to sale by or on the order of a licensed physician or other healthcare practitioner
	Product conforms to the Medical Device Directive 93/42/EEC (BSI)

Gamete

Oocyte Retrieval Buffer

REF ORB-20

REF ORB-50

Intended Use

Oocyte Retrieval Buffer is used in the removal of oocytes from ovarian follicles and to reduce stress on the oocytes during this procedure.

General Information

Oocyte Retrieval Buffer is supplemented with gentamicin (0.01 mg/mL).

pH range (25°C in air) is 7.30-7.50

Osmolality Range (mOsm/kg) is 280-290

Storage and Stability

Oocyte Retrieval Buffer must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Warm the Oocyte Retrieval Buffer to 37°C prior to use.

Sperm Buffer

REF SPB-20

REF SPB-50

Intended use

Sperm Buffer is used to wash and resuspend sperm for the insemination step in intrauterine insemination (IUI), IVF or in diagnostic washing.

General Information

Sperm Buffer is supplemented with human serum albumin (10 mg/mL) and gentamicin (0.01 mg/mL).

pH range (25°C in air) is 7.30-7.50

Osmolality Range (mOsm/kg) is 285-295

Storage and Stability

Sperm Buffer must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Warm the Sperm Buffer to 37°C prior to use.
- Sperm Buffer is suitable for both open and micro culture. If using micro drops, ensure culture oil is used to avoid evaporation and resultant osmotic stress.

Sperm Wash Gradient Set

REF SWG-01

Intended use

Sperm Wash Gradient Set is used to separate sperm from seminal plasma as well as separating highly motile sperm in preparation for insemination.

General Information

The Sperm Wash Gradient Set consists of 2 solutions:

- Sperm Wash Gradient (45%)
 - pH range (25°C in air) is 7.30-7.50
 - Osmolality Range (mOsm/kg) is 285-295
- Sperm Wash Gradient (90%)
 - pH range (25°C in air) is 7.30-7.50
 - Osmolality Range (mOsm/kg) is 280-290

The solutions are supplemented with human serum albumin (10mg/mL), gentamicin (0.01mg/mL) and coated silica.

Storage and Stability

Sperm Wash Gradient Set must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Equilibrate both 45% and 90% gradient solutions in the Sperm Wash Gradient Set to room temperature prior to use.
- Allow the semen to liquefy at 37°C for approximately 30 minutes.
- Prepare by adding 1.5 mL of Sperm Wash Gradient (90%) and then gradually overlaying 1.5 mL of Sperm Wash Gradient (45%) in conical base tubes.

- Layer up to 1.5 mL of the ejaculate on the gradient.
- 1st spin. Centrifuge 20 minutes at 300 g then carefully remove the seminal plasma, the upper interface, the 45% layer and the lower interface. Leave the remainder of the 90% layer intact.
- Remove the sperm pellet, place into a sterile conical tube and re-suspend in Sperm Buffer (SPB-20, SPB-50) (≥ 3 ml, ≤ 10 ml).
- 2nd spin. Centrifuge 6 minutes at 500 g. Remove the supernatant and re-suspend the pellet in a small volume (approximately 200 µL) of Sperm Buffer (SPB-20, SPB-50) or Sperm Medium (SPM-20, SPM-50).
- Assess the sperm preparation for motile concentration. Adjust as required.
- If sperm has been resuspended in Sperm Buffer, store in an incubator at 37°C until required.
- If sperm has been resuspended in Sperm Medium, store in a 6% CO₂ incubator at 37°C until required.

Growth

Sperm Medium

REF SPM-20

REF SPM-50

Intended Use

Sperm Medium is used to wash and resuspend sperm for the insemination step in IUI, IVF or in diagnostic washing, optimised for storage in a 6% CO₂ incubator.

General Information

Sperm Medium is supplemented with human serum albumin (10 mg/mL) and gentamicin (0.01 mg/mL).

pH range (37°C with 6% CO₂) is 7.30-7.50

Osmolality Range (mOsm/kg) is 285-295

Storage and Stability

Sperm Medium must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Semen should be allowed to liquefy at 37°C for approximately 30 minutes prior to use.
- Warm Sperm Medium at 37°C in a 6% CO₂ incubator prior to use.

For the swim-up technique:

- Gently underlay aliquots (300-1000 µL) of fully liquefied semen under 1.5 mL aliquots on equilibrated Sperm Medium in 14 mL conical tubes.
- Incline the tubes at an angle of 45°, to increase surface area of the semen-culture medium interface, and incubate in a 6% CO₂ incubator at 37°C for 60 minutes.
- Gently return tube to upright position remove the uppermost medium above the semen (approximately 0.7-1.0 mL).
- Assess the sperm preparation for motile concentration. Adjust as required.
- Store in a 6% CO₂ incubator at 37°C until required.

Gems Fertilisation (FEM-20, FEM-50), Cleavage (CLM-20, CLM-50) and Blastocyst (BLM-20, BLM-50) Media are optimised for use in a low oxygen (5%) environment.

Fertilisation Medium

REF FEM-20

REF FEM-50

Intended Use

Fertilisation Medium is used to provide a suitable environment for both oocytes and sperm, to promote optimal fertilisation rates.

General Information

Fertilisation Medium is supplemented with human serum albumin (5 mg/mL) and gentamicin (0.01 mg/mL).

pH range (37°C with 6% CO₂) is 7.30-7.50

Osmolality Range (mOsm/kg) is 295-305

Storage and Stability

Fertilisation Medium must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Fertilisation Medium is suitable for both open and micro culture. If using micro drops ensure culture oil is used to avoid evaporation and resultant osmotic stress.

- A suitable wash volume must also be prepared.
- Warm the Fertilisation Medium to 37°C in a 6% CO₂ incubator for a minimum of 4 hours prior to use.
- Fertilisation Medium is suitable for a standard sperm exposure (16-18 hours) or a short insemination.
- Following the oocyte and sperm preparation, add the appropriate number of sperm to each well or droplet containing the oocytes and return the dish to the incubator until fertilisation check.
- Subsequent zygotes should be transferred to Cleavage Medium (CLM-20, CLM-50).

- After achieving desired developmental stage on Day 3, transfer the embryo(s) to a culture dish containing the pre-warmed Blastocyst Medium. Wash the embryo before placing it in the final drop/well and returning the dish to the incubator.
- The reproductive scientist may continue to grow the embryo until transfer or cryopreservation on Day 4 to 6.

Cleavage Medium

REF CLM-20

REF CLM-50

Intended Use

Cleavage Medium is used to provide necessary nutrients for embryo development *in vitro*. The embryos can remain in this solution for up to 2 days prior to being transferred to the uterus or grown for up to another 3 days in Blastocyst Medium (BLM-20, BLM-50).

General Information

Cleavage Medium is supplemented with human serum albumin (5 mg/mL) and gentamicin (0.01 mg/mL).

pH range (37°C with 6% CO₂) is 7.30-7.50

Osmolality Range (mOsm/kg) is 285-295

Storage and Stability

Cleavage Medium must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Cleavage Medium is suitable for both open and micro culture. If using micro drops ensure culture oil is used to avoid evaporation and resultant osmotic stress.
- A suitable wash volume must be prepared to wash the embryo from the Fertilisation Medium (FEM-20, FEM-50).
- Warm Cleavage Medium in a 6% CO₂ incubator at 37°C for a minimum of 4 hours prior to use.
- Following fertilisation, transfer the embryos from the Fertilisation Medium to a culture dish containing the pre-warmed Cleavage Medium. Wash the embryo before placing it in the final drop/well and returning the dish to the incubator.
- The reproductive scientist may continue to grow the embryo until transfer on Day 3 or transfer to Blastocyst Medium (BLM-20, BLM-50) for further growth from Day 3 up to Day 6.

Blastocyst Medium

REF BLM-20

REF BLM-50

Intended Use

Blastocyst Medium is used to provide necessary nutrients for embryo development *in vitro* from the cleavage stage to the blastocyst stage. The embryos may remain in this solution for up to 3 days prior to being transferred to the uterus.

General Information

Blastocyst Medium is supplemented with human serum albumin (5 mg/mL) and gentamicin (0.01 mg/mL).

pH range (37°C with 6% CO₂) is 7.25-7.45

Osmolality Range (mOsm/kg) is 285-295

Storage and Stability

Blastocyst Medium must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Blastocyst Medium is suitable for both open and micro culture. If using micro drops ensure culture oil is used to avoid evaporation and resultant osmotic stress.
- A suitable wash volume must be prepared to wash the embryo from the Cleavage Medium (CLM-20, CLM-50).
- Warm Blastocyst Medium in a 6% CO₂ incubator at 37°C for a minimum of 4 hours prior to use.

Gerri Medium

REF ONE-20

REF ONE-50

Intended Use

Gerri Medium is used to provide necessary nutrients for embryo development *in vitro*. The embryos can remain in this solution for up to 5 days prior to being transferred to the uterus.

General Information

Gerri Medium is supplemented with human serum albumin (5 mg/mL) and gentamicin (0.01 mg/mL).

pH range (37°C with 6% CO₂) is 7.30-7.50

Osmolality Range (mOsm/kg) is 285-295

Storage and Stability

Gerri Medium must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used.
- Gerri Medium is suitable for both open and micro culture. If using micro drops ensure culture oil is used to avoid evaporation and resultant osmotic stress.
- A suitable wash volume must be prepared to wash the embryo from the Fertilisation Medium (FEM-20, FEM-50).
- Warm Gerri Medium in a 6% CO₂ incubator at 37°C for a minimum of 4 hours prior to use.
- Following fertilisation, transfer the embryos from the Fertilisation Medium to a culture dish containing the pre-warmed Gerri Medium. Wash the embryo before placing it in the final drop/well and returning the dish to the incubator.

Cryopreservation

Gavi Medium Cartridge

REF GAVI-MED-20

Intended Use

Gavi Medium Cartridge is intended for use within the Gavi instrument (GAVI-INS-01) for the vitrification of embryos for Assisted Reproductive Technology (ART) procedures.

General Information

The Gavi Medium Cartridge consists of 2 solutions:

- Gavi solution 1 which has a white twist-top cap is the equilibration solution supplemented with human serum albumin (16.8 mg/mL),
- Gavi Solution 2 which has a red twist-top cap is the vitrification solution supplemented with human serum albumin (13.5 mg/mL).

pH range (25°C in air) is 7.30-7.50

Gavi Medium Cartridge is for single use only.

Storage and Stability

Gavi Medium Cartridge must be stored in the original clear plastic packaging tray at all times, prior to use to ensure traceability and refrigerated at 2-8°C. Do not freeze.

Directions for Use

Refer to Gavi User Manual, section 5.2 Preparing for Consumables and Accessories, on use of the Gavi Medium Cartridge with the Gavi instrument. If the Gavi User Manual has been misplaced, please contact Genea Biomedx for a replacement at info@geneabiomedx.com

Vitrification Set

REF VIT-01

Intended Use

Vitrification Set is used for the vitrification of embryos for Assisted Reproductive Technology (ART) procedures.

General Information

The Vitrification Set consists of 3 solutions:

- Vitrification Set VitSol 1 – supplemented with human serum albumin (18.3 mg/mL)
- Vitrification Set VitSol 2 – supplemented with human serum albumin (16.2 mg/mL)
- Vitrification Set VitSol 3 (DMSO)

pH range (25°C in air) is 7.30-7.50

This Set is designed to complement the Gems Warming Set (WRM-01).

⚠ The Vitrification Set is not ready for use until the 3 solutions have been combined as per the Directions for Use. VitSol 3 is Dimethyl sulphoxide (DMSO). A Material Safety Data Sheet (MSDS) is available upon request. ⚠ The Vitrification Set must be used with a legally marketed storage device that is indicated for use in vitrification procedures. The storage device must be a closed system to prevent the potential risk of viral contamination.

VitSol 1 and VitSol 2 are supplemented with human serum albumin.

Storage and Stability

VitSol 1 and VitSol 2 in the Vitrification Set must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

⚠ VitSol 3 (DMSO) should be removed from the Vitrification Set packaging and stored at room temperature in original unopened container and protected from light.

Directions for Use

Preparation

- Aseptic technique should be used.
- DMSO freezes at <18.5°C. Equilibrate VitSol3 (DMSO) to room temperature prior to use.
- It is recommended that solutions with DMSO added do not sit for more than 1 hour. Discard after use.

All vitrification solutions should be prepared in an embryo safe vessel (for example in individual wells of a 4-well dish or in a small tube) just prior to use. Solutions should be prepared in the following ratios:

- For Vitrification Solution 1: Add 40 µL of VitSol 3 (DMSO) for every 460 µL of VitSol 1
- For Vitrification Solution 2: Add 80µL of VitSol 3 (DMSO) for every 420 µL of VitSol 2

Care should be taken when measuring and adding VitSol 3 (DMSO) into VitSol 1 & VitSol 2, as it is a viscous solution that requires careful measuring and thorough mixing, for example by pipetting the solution in and out a few times after addition.

Method

- Prepare the vitrification device for use as per manufacturer's instructions.
- Place the embryo to be vitrified into Vitrification Solution 1 for equilibration (collapse and re-expansion to 80% of original volume).
- When equilibration is complete, transfer the embryo to a well with a minimal amount of Vitrification Solution 2.
- Thoroughly wash the embryo. Transfer the embryo to the vitrification device and vitrify according to the Instructions for Use of the vitrification carrier. Time from transfer to Vitrification Solution 2 until vitrified should not exceed 90 seconds.
- Following vitrification, transfer the vitrified cells to the storage tank. Make sure the vitrified cells are submerged under liquid nitrogen at all times.

Warming Set

REF WRM-01

Intended Use

Warming Set and its individual solutions (WarmSol 1, WarmSol 2, WarmSol 3) is used for the warming of embryos that have undergone vitrification using the Vitrification Set (VIT-01) or Gavi Medium Cartridge (GAVI-MED-20). WarmSol 1, WarmSol 2 and WarmSol 3 should not be used individually, always use as part of the Warming Set.

General Information

The Warming Set consists of 3 solutions:

- Warming Set WarmSol 1 – supplemented with human serum albumin (20.0 mg/mL)
- Warming Set WarmSol 2 – supplemented with human serum albumin (20.0 mg/mL)
- Warming Set WarmSol 3 – supplemented with human serum albumin (20.0 mg/mL)

The solutions of the Warming Set are supplemented with human serum albumin.

pH range (25°C in air) is 7.30-7.50

Osmolality (mOsm/kg)

- Warming Set WarmSol 1 – 1280-1320
- Warming Set WarmSol 2 – 780-820
- Warming Set WarmSol 3 – 295-305

Storage and Stability

The Warming Set must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

Preparation

- Aseptic technique should be used.
- Equilibrate WarmSol 1 to 37°C before use.
- Equilibrate WarmSol 2 & WarmSol 3 to room temperature before use.
- For the culture of the embryo post warming, prepare a suitable volume of appropriate culture medium in a 6% CO₂ incubator at 37°C for a minimum of 4 hours prior to use.

Method

- Identify the embryo and vitrification device to be warmed as per the manufacturer's instructions.
- Add 2.0mL WarmSol 1 to a small petri dish and equilibrate to 37°C prior to use.
- Using a 4 well dish, add 500µL WarmSol 2 into well 1 and 500µL WarmSol 3 into wells 3 and 4. Equilibrate to room temperature.
- Extract the embryo from the vitrification device and place into the small petri dish with WarmSol 1.
- After 1 minute transfer the embryo from the petri dish to well 1 with WarmSol 2 for 3 minutes to begin re-equilibration.
- Transfer the embryo to well 3 with WarmSol 3 for 5 minutes.
- Transfer the embryo to well 4 with WarmSol 3 for 1 minute. Wash the embryo well.
- Transfer the embryo to a dish containing the appropriate equilibrated culture medium and incubate until transfer.

Other

VitBase

REF VBS-20

Intended Use

VitBase is a HEPES buffered medium that is used to maintain embryos for a short period of time in a non-gassed environment.

General Information

Can be used at either room temperature or after equilibration to 37°C.

VitBase is supplemented with human serum albumin (20 mg/mL).

pH range (25°C in air) is 7.30-7.50

Osmolality Range (mOsm/kg) is 295-305

Storage and Stability

VitBase must be stored in original unopened container and protected from light, refrigerated at 2-8°C. Do not freeze.

Directions for Use

- Aseptic technique should be used
- Aliquot the required amount into a suitable embryo safe vessel and either allow it to come to room temperature or equilibrate to 37°C before use.

MANUFACTURER

Genea Biomed Pty Ltd
Level 2, 321 Kent Street, Sydney
NSW, 2000, AUSTRALIA

Email: info@geneabiomedx.com

Web: www.geneabiomedx.com

AUTHORISED REPRESENTATIVE – EUROPE

EC REF

DONAWA LIFESCIENCE CONSULTING SRL

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