



Guardian

The electronic witnessing system using dual identification label technologies, designed for simplified sample traceability for patients and IVF clinics.



INTRODUCING GUARDIAN

Guardian is an electronic witnessing system using both types of identification label technologies, RFID chip and QR codes, designed as a simplified sample traceability system for patients and IVF clinics. The customer can choose the label options for each ID check procedure to achieve an optimized workflow. The hand-held PDA can read both types of labels which allows for ID checks to be done anywhere and everywhere. Guardian aims to help labs manage daily workflow and ensure ID checks are performed per protocol. Guardian, offering seamless traceability everywhere you need it.



Customizable System

It seamlessly adapts to the workflow of your lab, rather than requiring your lab to adapt to it.



Seamless Integration


Hardware has a small footprint, is portable, and enables seamless implementation.



Diverse Workflow Display

ID-check status and workflow can be viewed from day zero to day seven.

GUARDIAN: SEAMLESS TRACEABILITY, EVERYWHERE YOU NEED IT





Personal Digital Assistant (PDA)

Witnessing on the move for RFID and QR codes.



Workstation

Guardian can be used in desktop mode in a network connected environment.



RFID All-in-one

Touch screen workstation with RFID reading built-in*.



Monitor

Project guardian to an external lab display to quickly monitor progress and alarms.



Label Printer

Enabling printing of RFID tags, QR, and cryogenic labels for your workflow.



Server

The beating heart of guardian containing software & databases, with connectivity to the clinic network.

*Note: Final appearance of interfaces and devices may vary; subject to final release or customization.

GUARDIAN: VERIFIES PATIENT ID ACROSS ALL STEPS AND LOCATIONS



Seamless Traceability


Write electronic information and print visual labels simultaneously on RFID tags.

Flexibility Delivered

Select RFID or QR codes where needed.

Undisruptive Installation

No physical hardware change to existing equipment and lab set up.




Mobility Built-in

Mobile witnessing anytime and everywhere.



Transparency Connected

Mismatches distributed to PDA and dashboard for full visibility.



For healthcare professionals only. Please refer to the instructions for use.
For further information, please contact your Service Representative or visit:
www.geneabiomedx.com