



Guardian

The electronic witnessing system using dual identification label technologies, designed for seamless protection for patients and IVF clinics.



Patient List Batch Printing Comprehensive Enquiry Statistical Analysis Consumable Management System Configuration
Jason Taylor Embryologist

Active Cycle Full Process Mode Console Mode

Remaining 457 cycles

Incoming	Day 0	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
Cash Virginia 3232 Fresh Cycle, 3-11	Kingsley Bronwyn 1387 Fresh Cycle, 11-20	Ball Lucille 1985 Thaw Cycle, 3-13	Duck Daisy 5683 Fresh Cycle, 1-4	Ball Lucille 1985 Thaw Cycle, 3-5	Brown Corrine 0999 Thaw Cycle, 11-5	Carr Alicia 12366 Thaw Cycle, 2-3	Doe Jane 9285 Thaw Cycle, 4-3	Edwards Katrina 0962 Thaw Cycle, 10-1
McDonald Sydney 4863 Fresh Cycle, 1-12	Blackburn Jasna 2784 Fresh Cycle, 11-21	Anderson Apple 6314 Thaw Cycle, 3-12	Blackburn Jasna 2784 Fresh Cycle, 11-3	Albright Sally 1765 Thaw Cycle, 3-1	Clake Amelia 3542 Thaw Cycle, 2-2	Casir Belicia 1589 Thaw Cycle, 11-10	Estevez Maria 9866 Thaw Cycle, 10-2	
Duck Daisy 5683 Fresh Cycle, 1-23	Albright Sally 1765 Thaw Cycle, 3-14	Cook Lily A40261 Fresh Cycle, 20240391-1	Anderson Apple 6314 Thaw Cycle, 3-2					
Olsson Sandy 1236 Fresh Cycle, 11-22	Bening Becky 7887 Thaw Cycle, 3-11					Brown Melissa 2546 Thaw Cycle, 2-1		

INTRODUCING GUARDIAN

Guardian is an electronic witnessing system using both types of identification label technologies, RFID chip and QR codes, designed as a seamless protection system for patients and IVF clinics. The customer can choose the label options for each ID check procedure to achieve an optimized workflow. The hand-held PDA can read both types of labels which allows for ID checks to be done anywhere and everywhere. Guardian aims to help labs manage daily workflow and ensure ID checks are performed per protocol. Guardian is your reliable lab safety management assistant.



Customizable System

It seamlessly adapts to the workflow of your lab, rather than requiring your lab to adapt to it.



Seamless Integration


Hardware has a small footprint, is portable, and enables seamless implementation.



Diverse Workflow Display

ID-check status and workflow can be viewed from day zero to day seven.

GUARDIAN: SEAMLESS SAFETY, EVERYWHERE YOU NEED IT





Personal Digital Assistant (PDA)

Witnessing on the move for RFID and QR codes.




Workstation

Guardian can be used in desktop mode in a network connected environment.



RFID All-in-one

Touch screen workstation with RFID reading built-in.



Monitor

Project guardian to an external lab display to quickly monitor progress and alarms.



Label Printer

Enabling printing of RFID tags, QR, and cryogenic labels for your workflow.



Server

The beating heart of guardian containing software & databases, with connectivity to the clinic network.

GUARDIAN: VERIFIES PATIENT ID ACROSS ALL STEPS AND LOCATIONS



Advanced Safety


Write electronic information and print visual labels simultaneously on RFID tags.

Flexibility Delivered

Select RFID or QR codes where needed.

Undisruptive Installation

No physical hardware change to existing equipment and lab set up.




Mobility Built-in

Mobile witnessing anytime and everywhere.



Transparency Connected

Mismatches distributed to PDA and dashboard for full visibility.



For healthcare professionals only. Please refer to the instructions for use.
For further information, please contact your Service Representative or visit:
www.geneabiomedx.com